

BATERIAS DE USO MEDICO

ALARIS MEDICAL

Modelo	Voltaje / Capacidad
230, 530, 600 Volumetric Infusion Pump/Controller	12 Volt / .70 Ah
260, 260+, 260IR, 262, 262A, 560, 560+, 565, 570 Pump	12 Volt / 2.5 Ah
280, 580, 580 Plus, 590, 599 Star-Flow Volumetric Pump, 3080 Keofeed Enteral Feeding Pump	6 Volt / 3.4 Ah
4000, 4000A, 4000EE, 4200, 5200 Vital Check Monitor	12 Volt / 2.5 Ah
4400, 4410, 4415, 4510 Vital Check Monitor	6 Volt / 4.5 Ah
500, 3000 Keofeed Enteral Feeding Pump	12 Volt / 2.5 Ah
700, 710 Syringe Infusor	8 Volt / 2.5 Ah

BIO - MEDICAL DEVICES

Modelo	Voltaje / Capacidad
Crossvent 2i (PRT4402)	7.2 Volt / 3.0 Ah
Crossvent 2i+	7.2 Volt / 4.6 Ah
Crossvent 3+, 4+ (PRT4467) (NiMH - 3 wire)	7.2 Volt / 10.0 Ah
Crossvent 3, 3+, 4, 4+ (NiCd - 2 wire)	7.2 Volt / 5.0 Ah
M10 Ventilator Monitor (Save Original Wire Harness)	7.2 Volt / .85 Ah

CARDIOLINE

Modelo	Voltaje / Capacidad
AR1200 ECG	12 Volt / 1.6 Ah
AR2100 ECG	12 Volt / 2.5 Ah
AR600 ECG	9.6 Volt / 1.6 Ah
Delta 1 Plus, Delta 3 Plus (1220211-01)	24 Volt / .70 Ah
Delta 60 Plus	24 Volt / 1.8 Ah
EP700 EKG	18 Volt / .70 Ah
Recorder-E-1, E3 ECG (Recargable) (Requiere 8/unidad)	1.2 Volt / 2.3 Ah

CAREFUSION

Modelo	Voltaje / Capacidad
Aveo Ventilator 68339A	24 Volt / 4.5 Ah
LTV1150, LTV1200 Sprint Pack Battery (19444-001) (OEM)	14.4 Volt / 6.6 Ah
NOX-T3 Portable Sleep Monitor	1.2 Volt / 2.1 Ah

CRITICARE SYSTEMS, INC.

Modelo	Voltaje / Capacidad
508 EEG Monitor	12 Volt / 3.8 Ah
8100 Muliti Parameter (Requiere 2/unidad)	12 Volt / 2.3 Ah
8100E, 8100EP, 8100EP1 (Requiere 2/unidad)	12 Volt / 5.0 Ah
8100E1, nGenuity, nCompass vital signs monitor	12 Volt / 5.0 Ah
ELB External Power Supply (2035614-001)	24 Volt / 4.0 Ah
Oxygen Blood Pressure Unit 507N	12 Volt / 2.5 Ah
Poet I, Poet II, Poet IQ	12 Volt / 3.4 Ah
Poet IQ2, (8500Q) CO2 Monitor (42381B001)	12 Volt / 2.15 Ah
Poet LT Pulse Oximeter, 507 NJC BP, 602-14 (80512B001)	7.2 Volt / 2.7 Ah
Poet Plus 8100 (83278B001) *Short Version* 5.65" (Requiere 2/unidad) * May use either Short or Long Version*	12 Volt / 2.0 Ah
Poet Plus 8100 *Long Version* 7.17" (Requiere 2/unidad)	12 Volt / 2.3 Ah

Modelo	Voltaje / Capacidad
Poet TE Plus Oximeter, 602-11, 602-13	12 Volt / 2.5 Ah
Pulse Oximeter 500	12 Volt / 1.3 Ah
Pulse Oximeter 501	4.8 Volt / .70 Ah
Pulse Oximeter 502, 504, 504US, 50245, 50445, 506, DXN	6 Volt / 4.5 Ah
Pulse Oximeter 503, 503SP02T, 1040, 5040, Poni	6 Volt / 1.2 Ah
Pulse Oximeter 503S	6 Volt / 1.2 Ah
Pulse Oximeter 504DX	6 Volt / 3.2 Ah
Pulse Oximeter 506NT3, VSM 506 Series	6 Volt / 7.0 Ah
Pulse Oximeter 600, 601, 602, 602-1, 602-2, 602-3	12 Volt / 3.4 Ah
Pulse Oximeter END/TLCO2	12 Volt / 3.4 Ah
Vital Signs Monitor 503, 507, 1100	12 Volt / 3.4 Ah
Vital Signs Monitor 507S	12 Volt / 3.8 Ah
VSM 5070P, 507P, 507ER, 507N, 507E, Scholar II, III, EL	12 Volt / 2.5 Ah

DRAEGER

Modelo	Voltaje / Capacidad
Apollo (8603057) (Set de 2 baterías)	24 Volt / 17.0 Ah
Baby Log 8000	8.4 Volt / .180 Ah
Babylog VN500 Ventilator (8415290-08) (Send in for retrofit)	24 Volt / 3.0 Ah
CPU Board (Memory)	3 Volt / 1.5 Ah
Desflurane Vaporizer (M35526)	6 Volt / .51 Ah
Evita (Memory)	3.6 Volt / .17 Ah
Evita Dura Ventilator (Externa)	12 Volt / 33.0 Ah
Evita XL (Memory)	3 Volt / 1.5 Ah
Evita, Savina XL, 2, 4 Ventilator (Interna) (Panasonic) (2/unidad)	12 Volt / 3.4 Ah